

# FACT SHEET



**CS 7200**  
Imaging Plate System

# CS 7200 Imaging Plate System

## Product Summary

The CS 7200 is Carestream Dental's ultra-compact imaging plate system that delivers leading image quality. The slim design and the quiet scanning of the unit are ideal for chairside use and the system features a familiar film-like workflow, which paves the way for practices looking to transition to digital in a cost-effective manner.

## Product Features and User Benefits

**Exceptional Image Quality in Less Time:** With 17 line pairs/mm true resolution, the CS 7200 delivers images that are among the highest resolution available in just eight seconds for fast, confident diagnoses every time. New image processing delivers sharp, optimized contrast and low noise images. Pre-programmed anatomic filters, including endodontic, periodontic and dentin-enamel junction filters, enhance image contrast according to diagnostic needs with one click. Also, three scanning modes allow the user to give priority to either speed or resolution to best match diagnostic needs.

**Thin, Flexible Plates:** The CS 7200 features flexible and comfortable plates—available in sizes 0, 1 and 2—that are ideal for daily intraoral indications such as periapical, bitewing and pediatric exams. These durable, scratch-resistant plates require no chemicals and can be used hundreds of times, therefore, eliminating the cost of film and processing consumables. In fact, newly designed hygienic sheaths facilitate plate insertion while minimizing risk of cross-contamination. Afterwards, plates are ready for reuse immediately after processing.

**A Film-like Workflow with the Benefits of Digital:** Imaging plate technology provides many practitioners a better compromise between film workflow and digital benefits. With its similar handling to analog film, the CS 7200 offers a smoother transition to digital. The system is a great solution for renewing or transitioning to digital equipment, which means limited training and less of a learning curve when making the switch.

**User-Friendly Operation:** Operating the new CS 7200 is easy and user-friendly. With new direct acquisition, an image appears directly in the opened imaging software, with no need to handle any acquisition interface or manipulate a mouse. To avoid the risk of human error and lost images, the unit doesn't process a plate if the plate is oriented upside down.

**Cost-Effective and Compact:** The compact and space-saving design—just 13 cm wide—and quiet operation of the CS 7200 is perfect for chairside use or small practices. Not only that, its affordable price, reliability and ease of maintenance make it a cost-effective digital imaging solution with no compromise on quality.



**Minimal Training Required:** The CS 7200 requires minimal training or disruption to established work patterns. The unit connects directly to the computer with a standard USB connection. Carestream Dental's imaging software can be easily integrated with leading dental practice management software, thereby minimizing the learning curve for your staff. Plus, many services (configuration, diagnosis, repair, calibration, software updates) can be performed remotely, reducing downtime and keeping service cost under control.

### More Information

For more information on Carestream Dental's innovative solutions or to request a product demonstration, call **800.944.6365** or visit **www.carestreamdental.com**.

### About Carestream Dental


Carestream Dental provides industry-leading imaging, CAD/CAM, software and practice management solutions for dental and oral health professionals. With more than 100 years of industry experience, Carestream Dental products are used by seven out of 10 practitioners globally and deliver more precise diagnoses, improved workflows and superior patient care. For more information or to contact a Carestream Dental representative, call **800.944.6365** or visit **www.carestreamdental.com**.

### About Carestream Health

Carestream is a worldwide provider of dental and medical imaging systems and IT solutions; X-ray imaging systems for non-destructive testing and advanced materials for the precision films and electronics markets. For more information about the company's broad portfolio of products, solutions and services, please contact your Carestream representative, call **888.777.2072** or visit **www.carestream.com**.

**PR Contact:** Amanda Gong **770.226.3213**, [amanda.gong@carestream.com](mailto:amanda.gong@carestream.com)

### Follow Carestream Dental online



©Carestream Health, Inc. 2016. 14037 DE CS 7200 0516