

FACT SHEET



CS 1200
Intraoral Camera

CS 1200 Intraoral Camera

Product Summary

The CS 1200 intraoral camera serves as an affordable entry point into digital imaging for almost any practice. In addition to Carestream Dental's industry-leading image quality, the CS 1200 is easy to share and integrate into existing workflows.

Product Features and User Benefits

Industry-Leading Image Quality: The CS 1200 offers a still image resolution of 1024 x 768 pixels for maximum detail and a video image resolution of 640 x 480 pixels. The camera's fixed focus range (3mm to 25mm) captures a variety of images. Additionally, the CS 1200 features a 6-LED auto-adjusting illumination system, which ensures perfectly and uniformly lit images.

Improved Patient Communication: The CS 1200 provides practitioners and patients with a detailed view of the patient's teeth, anatomical structures and mouth. Providing true-to-life images ultimately results in improved patient communication, increased case acceptance and better case documentation for sharing with third-parties and referrals.

Compact and Comfortable: The compact camera is lightweight and ergonomically designed to fit comfortably in hands of all sizes. With a conveniently placed acquisition button, the camera is easy to use, and practitioners can capture an image with one click. Built with a rounded head and tapered shape, the CS 1200 offers improved patient comfort during exams.

Convenient and Easy to Share: Designed with convenience and flexibility in mind, the CS 1200 offers USB and optional S-Video and AV connections to support both PC and analog displays. The camera also works without docking stations, making integration and sharing between operatories fast and convenient.

Built-In Memory: The CS 1200 can store up to 300 images within the camera itself. With a touch of a button, the images can be viewed on a TV, monitor or PC, eliminating the need for memory cards and computers in each operatory.

Affordable Entry Point into Digital Imaging: The CS 1200 is a high-quality intraoral camera at an affordable price-point that puts digital imaging within reach of almost any practice. The camera also comes with everything a practice needs to get started, including a mountable holder, hygienic sheaths and a USB cable.

Minimal Training Required: The CS 1200 is easy to learn and use, minimizing training time and improving practice productivity.



More Information

For more information on Carestream Dental's innovative solutions or to request a product demonstration, call **800.944.6365** or visit **www.carestreamdental.com**.

About Carestream Dental

Carestream Dental provides industry-leading imaging, software and practice management solutions for dental and oral health professionals. With more than 100 years of industry experience, Carestream Dental products are used by seven out of 10 practitioners globally and deliver more precise diagnoses, improved workflows and superior patient care. For more information or to contact a Carestream Dental representative, call **800.944.6365** or visit **www.carestreamdental.com**.

About Carestream Health

Carestream Health is a worldwide provider of dental and medical imaging systems and healthcare IT solutions; X-ray film and digital X-ray systems for non-destructive testing; and advanced materials for the precision films and electronics markets.

For more information about the company's broad portfolio of products, solutions and services, please contact your Carestream Health representative or visit **www.carestream.com**.

PR Contact: Amanda Gong **770.226.3213**, amanda.gong@carestream.com

Follow Carestream Dental online



©Carestream Health, Inc. 2016. 14035 DE 12 0516