

FACT SHEET



CS 7600
Imaging Plate System

CS 7600 Imaging Plate System

Product Summary

The CS 7600 is Carestream Dental's digital intraoral radiography system that reinvents imaging plate technology—improving usability, productivity and security with an innovative and fully automated workflow.

Product Features and User Benefits

Scan & Go Technology: The CS 7600 includes exclusive Scan & Go technology¹ that uses Smart Plates to prevent plate mix-up and reduce operation time. This feature—an industry first—identifies the plates electronically with patient information, allowing the device to automatically route scanned images to the proper computer and patient file.

Smart Plates: The CS 7600 uses thin, flexible Smart Plates that have no attached cable. The plates are available in size 0 to 4, and they are comfortable for patients and as easy to place as film.

User-Friendly Operation: The system's user-friendly operation automatically detects improper plate insertion, providing clear instructions for correction on a large-format, color display. The CS 7600 also supports batch scanning of mixed patient plates effortlessly and securely, removing the risk of operator mistakes. Practitioners can easily scan and display single or multiple plates, as well as full mouth series. Afterwards, plates are automatically erased for reuse. Multiple users can work on the system at the same time, as the system is always ready for operation.

Exceptional Image Quality—Fast: With up to 17 line pairs/mm true resolution, the CS 7600 delivers images that are among the highest resolution available for confident diagnoses every time. The unit is also among the fastest plate systems, with an image viewable in as little as five seconds.

Cost-Effective and Compact: Designed for chairside and shared use, the CS 7600 does not require a dedicated computer to run, which optimizes practices' investment and space. The unit's small footprint enables easy installation, wherever it will best serve patients—on a desktop or mounted to a wall.

Built-In Memory: The CS 7600 offers a built-in memory function to ensure images are captured and stored, even if there's a network failure during the process. The Scan & Go option automatically restores and sends the images to the appropriate computer when the network is repaired.

Minimal Training Required: As easy to use as film with all the benefits of digital imaging, the CS 7600 requires minimal training or disruption to established work patterns. Its images can be previewed on the unit, so users instantly know whether the desired image was successfully captured.



¹Available as optional feature on CS 7600.

More Information

For more information on Carestream Dental's innovative solutions or to request a product demonstration, call **800.944.6365** or visit **www.carestreamdental.com**.

About Carestream Dental

Carestream Dental provides industry-leading imaging, software and practice management solutions for dental and oral health professionals. With more than 100 years of industry experience, Carestream Dental products are used by seven out of 10 practitioners globally and deliver more precise diagnoses, improved workflows and superior patient care. For more information or to contact a Carestream Dental representative, call **800.944.6365** or visit **www.carestreamdental.com**.

About Carestream Health

Carestream Health is a worldwide provider of dental and medical imaging systems and healthcare IT solutions; X-ray film and digital X-ray systems for non-destructive testing; and advanced materials for the precision films and electronics markets.

For more information about the company's broad portfolio of products, solutions and services, please contact your Carestream Health representative or visit **www.carestream.com**.

PR Contact: Amanda Gong **770.226.3213**, amanda.gong@carestream.com

Follow Carestream Dental online



©Carestream Health, Inc. 2016. 14038 DE CS 7600 0516